

Genuine, washable

SHEEPSKIN & WOOL PILE



for
Medical/Home Health Care Use

Manufactured/distributed by:

The Wool Table Sheepskins

Workshop: 80 Barbados Blvd, Unit 9, Scarborough, ON M1J 1K9

Mailing: P.O. Box 88036, Scarborough, ON M1M 3W1

Tel: (416) 755-9103 (800) 854-9774

Fax: 647 351-5262

E-mail: thewooltable@rogers.com

Manufacturer/distributor of quality sheepskin products for more than 25 years.



MEDICAL SHEEPSKIN ACCESSORIES

*The Wool Table Sheepskins is a Canadian business manufacturing and distributing premium quality medical grade home health care products that will contribute to the comfort and relief of your customers. Product design is important to us. A major benefit of our products is **pressure relief** to reduce the risk of bed sores. Our sheepskin items are **machine washable** with the sheared wool length approximately 1-1.25" for optimum comfort and relief. The wool fibres permit the air to flow through them, **absorbing moisture** that is wicked to the atmosphere, providing a natural protection that is difficult to duplicate. Sheepskin is beneficial through **all** seasons. Please remember that sheepskin is a natural material that may contain minor irregularities or flaws. We do our best to remove them for you. The Wool Table Sheepskins is able to offer you & your customers three levels of wool materials – **wool pile, regular medical grade, and Hi-temp medical grade**. Thank you for your interest in our quality products that you have recognized for their natural ability to provide comfort and relief to convalescing people through **all** seasons. The merits of these products sell themselves. Custom work is available.*

Wool fleece colour code for sheepskins:

02 – cream, 07 – charcoal grey, 08 - black, 09 - pale gold.

WOOL & WOOL PILE

(Wool density more than 1150gms per sq. mtr)

N.B. Our wool pile material is "Down Under" Superwash 1" fleece with good density bonded into a polyester knit backing. Machine washable. Material made in USA. Colour: Natural.

Product No.	Description
1210	Bed pad , 20" x 24", soft Superwash 100% wool fleece pad.
1212	Bed pad, 24" x 30", as for #1210.
1214	Bed pad, 30" x 40", as for #1210.
1216	Bed pad, 30" x 60", as for #1210
1218	Bed pad, 30" x 72", as for #1210.
1221	Mattress pad, twin , 39" x 72", a plush 1" wool pile pad, c/w bias tape edging and elastics to hold @ the corners plus one centre strap.
1222	Mattress pad, double , 54" x 72", same as for #1221.
1223	Mattress pad, queen , 60" x 80", same as for #1221.
1224	Mattress pad, king , 78" x 72", same as for #1221.
1225	Mattress pad, hospital , 36" x 80", same as for #1221.
1226	Mattress pad, ¾ , 48" x 72", same as for #1221.
1227	Mattress pad, crib , 28" x 50", same as for #1221.
1230	Seat pad , 14w"x10", c/w lined underside, bias tape trim, & elastic loops to enable rolling it up & carrying it with you. Use for arena, bingo, car seat, church pew, etc.
1231	Seat pad, 16W" x 14", same as for #1230.
1233	Seat pad, 18W" x 16", same as for #1230.
1235	Seat pad, 18W" x 18", same as for #1230.
1236	Seat pad, 20"W x 18", same as for #1230.
1239	Walker handgrips pads, set of 2 pcs, c/w velcro tabs to hold.
1241	Wheelchair armrest pads , 11", set of 2, c/w velcro tabs to hold.
1242	Wheelchair armrest pads, 14", set of 2, c/w velcro tabs to hold.
1243	Wheelchair seat pad , 20" x 20", c/w 2 straps & velcro tabs to close under the seat. Bias tape trimmed.

- 1244 Wheelchair seat pad, 22" x 22", same as for #1243.
- 1245 Wheelchair seat pad, 16" x 16", same as for #1243.
- 1246 Wheelchair seat pad, 18"W x 16", same as for #1243.
- 1247 Wheelchair seat pad, 18" x 18", same as for #1243.
- 1248 Wheelchair seatback pad, 18"W x 17", c/w 2 velcroed loops to hang off the handles.
- 1249 Wheelchair **seat/seatback pad** (one piece), c/w attachments to hold (2 velcro loops to hold on to the handles & 3 sets of straps to hold elsewhere). Bias tape trimmed. Size: 18"W x 36"L.
- 1251 **Knee pads**, small, approx. 10"H x 8"W, side 'V' cut out, use in bed, while in a chair or wheelchair, velcro straps to attach, etc.
- 1253 Knee pads, medium, approx. 12"H x 10"W, same as for #1251.
- 1255 Knee pads, large, approx. 14"H X 12"W, same as for #1251.
- 1263 **Cushion cover**, 20" x 20", replacement cover w/black nylon skirt.
- 1264 **Protective sleeve**, wool fleece lined sleeve with cuffs, great to protect the mid arm or shins from injury, enabling healing.
 Sizes: S (9.5"L x 13.5" cir.), M (11"L x 15" cir.), L (12.5L" x 16.5" cir.), XL (12.5"L x 18" cir.).
- 1265 **Cushion cover**, 16" x 16", replacement cover w/black nylon skirt.
- 1266 Cushion cover, 18"W x 16", replacement cover w/black nylon skirt.
- 1267 Cushion cover, 18" x 18", replacement cover w/black nylon skirt.
- 1271 **Elbow** protectors, "bean shaped", comfortable 1" wool pile, c/w velcro straps to hold and no seam @ point of elbow.
- 1274 **Bedrest** pad, 12"Wx36"L, bias tape bound edges, same as for #1732.
- 1275 **Convalescent boots**, closed toes, above ankle cut, velcro tabs to close across the instep, sedentary use, non-skid sole. S, M, L, XL
- 1276 **Heel** protectors, comforting wool pile, c/w 2 velcro straps to hold and no seam across the heel. Same style as #1726.
- 1277 **Convalescent boots**, closed toes, above ankle cut, velcro tabs across the in-step, for sedentary use, pair. Sizes: S, M, L, XL.
- 1278 **Heel protectors**, extended forward, same as for #1728. M or L
- 1279 **Convalescent boots**, open toes, above ankle cut, velcro closure over the in-step, for sedentary use, pair. Sizes: S, M, L, XL.
- 1282 **Wool**, natural wool roping for packing between toes etc. 10gm bag.
- 1284 Wool, as for #1282. 30gm bag.
- 1290 **Recliner chair** seatcover, approx. 20"W x 52"L, c/w straps to hold.
- 1292 **Forearm wrap**, small, prevent the forearm sweating against another surface, 8"L x 10"W tapering to 6.5", c/w 3 velcro straps to hold.
- 1293 Forearm wrap, medium, prevent the forearm sweating against another surface, 9"L x 11"W tapering to 8", c/w 3 velcro straps to hold.
- 1294 Forearm wrap, large, prevent the forearm sweating against another surface, 10.5"L x 12"W tapering to 8", c/w 3 velcro straps to hold.
- 1296 **Hot water bottle cover**, the fleecy exterior cover allows heat to dissipate more evenly. C/w velcro tabs to close. Approx. 8"W x 10"L.
- 7130 **Woolskin** sheepskin detergent, a specialty Tea Tree oil based wash for sheepskin and wool with its anti-bacterial qualities plus its anti-dust mite abilities. 500 ml btl. Imported from Australia.
- 7422 **Socks**, a unique, luxury blend of possum fur/merino wool, & mulberry silk, these rib knit style socks are very warm and comfortable. Imported from New Zealand. Colours: Mocha, Brown, Black.
 Sizes: S (4-6), M (7-9), L (10-12), XL (13-14). Made in New Zealand.
- 7452 **Gloves**, a unique luxury blend of possum fur/merino wool, & mulberry silk, these light weight gloves will surprise you with their warmth. Colours: Mocha or Charcoal grey. Sizes: S, M, L. Made in New Zealand.

SHEEPSKIN PADS, regular grade

Pressure relief bed pads

- 1012 Lambskin, Infantcare, natural shape, wool sheared to 1", approx. 22" x 35". Wool density suitable for use with frailer adults too. Col: Maize.
These medical grade skins can be washed up to 40°C.
- 1116 Sheepskin, natural shape, approx. 26" x 41". Col: 09 only
- 1124 Sheepskin, extended, natural shape, approx. 26" x 54". Col: 09 only
- 1125 Sheepskin, extended, natural shape, approx. 26" x 60". Col: 09 only
- 1126 Sheepskin, extended, natural shape, approx. 26" x 72". Col: 09 only
- 2110 Lambskin, natural shape, approx. 8 sq ft, drape over the couch, car seat, wheelchair, etc. Cols: 02, 07, 08, 09.
The following sheepskin bed pads are custom made to cover the whole mattress area. Elastic corner straps plus a centre strap are attached. Medical grade. Col: 09 only.
- 1141 **Mattress overlay, single** size, 39" x 72".
- 1142 Mattress overlay, **double** size, 54" x 72".
- 1143 Mattress overlay, **queen** size, 60" x 80".
- 1144 Mattress overlay, **king** size, 78" x 72".
- 1145 Mattress overlay, **hospital** size, 36" x 80".
- 1146 Mattress overlay, **twin (3/4)** size, 48" x 72".
- 1147 Mattress overlay, **crib** size, 28" x 52".
- 1151 Bed pad, 20" x 30", custom shaped.
- 1152 Bed pad, 24" x 30", custom shaped.
- 1153 Bed pad, 30" x 36", custom shaped.
- 1154 Bed pad, 30" x 40", custom shaped.
- 1155 Bed pad, 30" x 48", custom shaped.
- 1156 Bed pad, 30" x 60", custom shaped.
- 1157 Bed pad, 30" x 72", custom shaped.
- 1158 Bed pad, 36" x 60", custom shaped.

Wheelchair Accessories

All colours available.

NB. *These pads are made to fit directly to the chair minus any additional cushions or padding. To cover this extra padding, let us know so we can provide additional strapping.*

- 1511 Armrest pad, approx. 4"W x 11"L, c/w velcro straps to attach.
- 1514 Armrest pad, approx. 4"W x 14"L, c/w velcro straps to attach.
- 1515 Seat pad, 16" x 16", c/w straps & velcro tabs to close under the seat.
- 1516 Seat pad, 18"W x 16", same as for #1515.
- 1517 Seat pad, 18" x 18", same as for #1515.
- 1518 Seat back pad, 18"W x 17", c/w velcroed loop attachments to go over the handle bars to hold. Hangs in place.
- 1519 HUGO transport chair seat back pad, 19"W x 17.5"H.
- 1520 Seat pad, 20" x 20", same as for #1515.
- 1522 Seat pad, 22" x 22", same as for #1515.
- 1523 HUGO transport chair seat pad, 19.5"W x 15.5"D.
- 1524 HUGO Navigator seat pad, 14.5"W x 13.5"D.
- 1526 Seat pad, 20"W x 18"D, same as for #1515.
- 1527 Seat/seatback pad, 18" x 36", c/w straps to hold.
- 1528 Hanger wraps, protects legs from banging on the hangers, 10"H.
- 1529 Calf protector, 15.5" x 8", w/straps to hold.
- 1532 Foot plate pads, 6" x 8", hook velcro strips attached, pair.
- 1533 Foot plate pads, 7" x 10", hook velcro strips attached, pair.
- 1534 Foot plate pad, 7.5" x 15", velcro straps attached.

Walker chair (Rollator) pads,

Comforting pads for the seat and seat back. The seat pads c/w straps & velcro

tabs to hold under the seat.

- 1542 Seat pad, 15" x 11", fits Nexus 3 et al, straps to hold,
1544 Seat pad, 17"W x 12.5", fits Dana Douglas 4200WX, straps to hold.
1546 Seat pad, 19"W x 13", fits Dana Douglas 4200SWX, straps to hold.
1555 Walker chair back strap pad, 1.75"W x 18"L, c/w 3 velcro tabs to hold.
Fits Opal, Riva, Legacy, etc. (Provides padding on the front side only).
1558 Walker chair rear bar pad, 5.5" x 11" cylindrical wrap, velcro to close.
1559 Rollator seatback pad, approx. 14.75" x 4.75", fits Airgo Comfort Plus,
Free Spirit Ergo Plus, etc
1561 Seat pad, 13.25" x 7.5", fits Legacy, etc. C/w straps & velcro to hold.
1562 HUGO Rolling walker seat pad (Standard), 12"x 12".
1563 HUGO Rolling walker seat pad, Elite, 14"W x 12"D.
1564 Seat pad, 15.50"W x 10", fits Evolution Horizon.
1565 Seat pad, 15.25"W x 9.75", fits Evolution Piper.
1567 Seat pad, 17.25"W x 9.50", fits Dana Douglas Infiniti 500.
1568 Seat pad, 14" x 14.25", fits Dana Douglas 4200.
1569 Seat pad, 11.75" x 11.75", fits Ergo, Free Spirit standard seat.
1571 **Cushion cover**, 16" x 16", replacement cover w/black nylon skirt.
Fits up to 3"D cushion.
1572 Cushion cover, 18"W x 16", replacement cover, same as for #1571.
1573 Cushion cover, 18" x 18", replacement cover, same as for #1571.
1574 Cushion cover, 20" x 20", replacement cover, same as for #1571
1581 **Recliner** chair seat cover, approx. 19"W x 52"L, c/w 2" wide elastic
straps to secure it to the upper and lower seat back.
1583 Recliner chair armrest covers, 13.5"W x 15"L, pair.
1585 Recliner chair armrest covers, 16"W x 18"L, pair.
1587 Recliner chair footrest pad, 13"D x 20"W, c/w straps to hold.
1591 **Headrest** lateral support cover, Otto Bock 430 N1=1, or N1=3.
1593 Headrest lateral support cover, Otto Bock H1=1, approx. 13" x 6".
1594 Headrest lateral support cover, Otto Bock C1=1, approx. 11" x 5.5".
1595 Headrest lateral support cover, Otto Bock C1=3, approx. 8" x 5".

Protective Aids

- 1715 Forearm **crutch** pads, set of 2 upper pads, fits to supports.
1716 Crutch pads, standard style, 2 handgrip pads only, c/w velcro closures.
1717 Crutch pads, standard, set of 4 pieces (2 x upper & 2 x handgrip pads),
c/w velcro closures.
1718 Crutch pads, standard style, 2 upper pads only, c/w velcro closures.
1719 **Walker** handgrip pads, approx. 4.75"W x 5.25"L, c/w velcro closures to
enable wrapping around handles. Fits quad canes, walkers, etc. Pair.
1719-1 Quad cane handgrip pad, same description as for #1719. One only.
1720 Walker handgrip pads, approx. 4.75"W x 5.50"L, same as for #1719.
1721 **Elbow** protectors, dense fleece for comfort, c/w 2 padded velcro straps
for closure. Main panel (under point of elbow) 3.5" x 8". Col: 09.
No seam @ point of elbow for greater comfort & relief.
1722 **Cane** handgrip pad, T-shaped cane handle, elastics to hold in place.
1723 **Tilt chair** handgrip pads, 4.75"W x 3.5"L, velcro closures, pair.
1726 **Heel** protectors, dense fleece for comfort, c/w 2 padded velcro straps
for closure over the instep. No heel seam. Approx. 5" forward. Col: 09.
1727 Heel protectors, extended above ankle (2"), c/w 3 padded straps.
1728 Heel protectors, extended forward under sole & sides beyond the
original model (#1726). Sizes: M (+3.5" forward), L (+5" forward), or
XL (+7" forward). Includes padded velcro straps to close.

1732 **Bedrest** pad, 12" x 36", a protective pad to lay across the bed
providing comfort & relief for under the heels, the buttocks, or

shoulders, etc. Has extended fabric tabs to tuck under the mattress.
Cols: 02, 07, 09.

- 1735 **Pavlik** harness wraps, set of 5 pcs – 2-1"x5.5", 2-2"x3.5", & 1-2"x3.9". Cols: 01 (white), 02, 09.
- 1741 **Knee** pads, small, approx. 10"H x 8"W, side 'v' cut out, c/w 2 velcro straps, use in bed, while in a chair or wheelchair, etc. pair. Col: 09.
- 1743 Knee pads, medium, approx. 12"H x 10"W, pair, same as for #1741.
- 1745 Knee pads, large, approx. 14"H x 12"W, pair, same as for #1741.
- 1751 **Forearm wrap**, small, a protective wrap 8"L x 10"W tapering to 6.5", c/w 3 velcro straps to hold. Helps prevent skin breakdown. Col: 09.
- 1752 Forearm wrap, medium, a protective wrap, 9"L x 11"W tapering to 8", c/w 3 velcro straps to hold. Helps prevent skin breakdown. Col: 09
- 1753 Forearm wrap, large, a protective wrap 10.5"L x 12"W tapering to 8", c/w 3 velcro straps to hold. Helps prevent skin breakdown. Col: 09.
- 1767 **Coccyx pads**, 2 – 6" x 6" dense medical grade sheepskin pads to provide pressure relief in the coccyx area. Col: 09.
- 1769 **Ear** cushion, 1.5" deep foam cushion w/sheepskin padding on top & fabric skirt, overall size 12" square, c/w center hole to allow for the ear's pressure relief. Wool colour: 09.
- 1782 **Palm pads**, 2 shaped pads to wrap through the palms to protect against encroaching arthritic fingers, etc. A velcro strap across the wrist holds each in place. One left and one right. Helps sweat issue. Col: 09.
- 1784 **Palm cushion**, double layer of sheepskin to help exercise the fingers. Use in the palm. An elastic band goes around the wrist.

Footwear

Our open toe slippers have an adjustable flap across the in-step, and are cut ankle high at the rear. This provides an option for those who don't want a fully enclosed slipper. Exterior colour is tan.

- 1838 **Slippers, open toes, ladies'**, c/w an adjustable flap across the in-step. Composition sole. Indoor or outdoor use. Sizes: 5, 6, 7, 8, 9, 10
- 1839 Slippers, open toes, men's, c/w an adjustable flap across the in-step. Composition sole. Indoor or outdoor use. Sizes: 7, 8, 9, 10, 11, 12, 13.

These slippers come with a sheepskin upper and either a cowhide leather sole or an indoor/outdoor composition sole. These slippers have been opened part way down the in-step to enable those with high in-steps, swollen feet, diabetics, etc, to enjoy this comfort. The rolled cuff around the ankle can be turned up for added warmth. These slippers have 2 velcro straps to close over the in-step and around the ankle. Intended for ambulatory uses. They can be rinsed or washed out if needed, using warm water & Woolskin or a mild detergent. The leather colour is tan.

- 1846 Slippers, ladies', **household style**, c/w indoor/outdoor composition sole, wider opening over the instep with velcro tabs to close. Sizes: 5, 6, 7, 8, 9, 10.
- 1847 Slippers, ladies', household style, c/w cowhide leather sole, wider opening over the instep, with velcro tabs to close. Sizes: 5, 6, 7, 8, 9, 10.
- 1848 Slippers, men's, household style, c/w hard indoor/outdoor sole, wider opening over the instep with velcro tabs. Sizes: 7, 8, 9, 10, 11, 12, 13. Sizes 14 & 15 by special order and add \$10.00).
- 1849 Slippers, men's, household style, c/w cowhide leather sole, wider opening over the instep with velcro tabs. Sizes: 7, 8, 9, 10, 11, 12, 13. (14 & 15 by special order and add \$10.00).

Convalescent boots: *Above the ankle cut, wider opening, c/w velcro tabs to close across the in-step and around the ankle. Intended for sedentary indoor use. The exterior leather is off-white/tan. Seam free across the rear of the heel.*

Sizes: Ladies – S (5-6), M (7-8), L (9-10).

Men's – S (7-8), M (9-10), L (11-12), XL (13-14), add \$5.00.

- 1851 **Convalescent boots**, ladies, open toes, pair. Sizes: S, M, L,
 1852 Convalescent boots, men's, open toes, pair. Sizes: S, M, L, XL.
 1855 Convalescent boots, ladies. closed toes, pair.
 1856 Convalescent boots, men's, closed toes, pair.
 1861 Convalescent boots, Hi-split style, closed toes, Ladies. Sizes: S, M, L.
 1862 Convalescent boots, Hi-split style, closed toes, Men's. Sizes: S, M, L, XL.
 1885 **Insoles**, ladies', lambskin fleece with leather backing. Wool length approx. 17-20mm. Col: white. Sizes: 5/6, 7/8, 9/10.
 1886 Insoles, men's, as for ladies. Sizes: 7/8, 9/10, 11/12, 13-14 (+\$2.00).
 1889 **Pressure pads**, 2 – 2" x 3" pads to alleviate pressure points in footwear e.g. on the in-step, on the ankle, under the heel, etc.
 3331 **Seatbelt pad**, 11" long, fits 2" wide strap, c/w velcro tabs to close.

Hi-Temp GREEN Medical Sheepskin

The Wool Table Sheepskins includes in its product line up this specialty grade of medical sheepskin. The Hi-temp skins can be washed up to 80° and the regular grade medical skins can be washed @ 40°. While beneficial for use in **institutional** settings, it can be used in home settings too. These skins can be sanitized and re-used by another client.

The product designs are the same as the regular medical ones, just a different material and pricing structure. See our main catalogue for those images

- 8118 **Sheepskin**, natural shape, approx. 41"L.,
 8124 **Bed pad**, natural shape, approx. 54"L
 8125 Bed pad, natural shape, approx. 60"L.
 8126 Bed pad, natural shape, approx. 72"L.
 8141 **Mattress pad**, single, 39" x 72".
 8142 Mattress pad, double, 54" x 72".
 8143 Mattress pad, queen, 60" x 80".
 8144 Mattress pad, king, 72" x 78".
 8145 Mattress pad, hospital, 36" x 80".
 8146 Mattress pad, twin (3/4), 48" x 72".
 8150 **Bed pad, custom**, 20" x 24"
 8151 Bed pad, custom, 20" X 30".
 8152 Bed pad, custom, 24" x 30".
 8153 Bed pad, custom, 30" x 36".
 8154 Bed pad, custom, 30" x 40"
 8155 Bed pad, custom, 30" x 48"
 8156 Bed pad, custom, 30" x 60"
 8157 Bed pad, custom, 30" x 72"
 8158 Bed pad, custom, 26" x 72".
 8159 Bed pad, custom, 36" x 60".
 8511 **Wheelchair armrest pads**, 4"W x 11"L, velcro straps, pair.
 8514 Wheelchair armrest pads, 4"W x 14"L, velcro straps, pair.
 8515 Wheelchair **seat pad**, 16" x 16", straps attach under seat.
 8516 Wheelchair seat pad, 18"W x 16", straps attach under seat.
 8517 Wheelchair seat pad, 18" x 18", straps attach under seat.
 8518 Wheelchair seatback pad, 18"W x 17", loops attach to handlebars.
 8520 Wheelchair seat pad, 20" x 20", straps attach under seat.
 8526 Wheelchair seat pad, 20"W x 18", straps attach under the seat.
 8527 Wheelchair seat/seat back pad, c/w straps to attach, 18" x 36".
 8531 **Knee pads**, small, approx. 10"H x 8"W, same as for #1741.
 8533 Knee pads, medium, approx. 12"H x 10"W, same as for #1741.
 8535 Knee pads, large, approx. 14"H x 12"W, same as for #1741.
 8541 Broda Tilt Recliner #785, armrest pads, 23" L.

- 8571 **Cushion cover**, 16" x 16", replacement cover w/black nylon skirt.
- 8572 Cushion cover, 18"W x 16", replacement cover w/black nylon skirt.
- 8573 Cushion cover, 18" x 18", replacement cover w/black nylon skirt.
- 8574 Cushion, 20" X 20", replacement cover w/black nylon skirt.
- 8716 **Crutch pads**, lowers only, pair.
- 8717 Crutch pads, 2 upper & 2 lower. (Set of 4 pcs.)
- 8718 Crutch pads, uppers only, pair.
- 8719 **Walker** handgrip pads, same as for #1719, pair.
- 8719-1 Quad cane handgrip pad, same as for #1719, one only.
- 8721 **Elbow protectors**, pair, same as for #1721.
- 8726 **Heel protectors**, pair, same as for #1726.
- 8727 Heel protectors, extended up, 2" higher, one size, pair.
- 8728 Heel protectors, extended forward, pair, Sizes: M, L, or XL.
- 8732 **Bedrest pad**, 12" x 36", same description as for #1732.
- 8751 **Convalescent boots**, ladies', open toes. Pair. S, M, L.
- 8752 Convalescent boots, men's, open toes, pair. S, M, L, XL. (add \$5.00).
- 8755 Convalescent boots, ladies', closed toes, pair. S, M, L.
- 8756 Convalescent boots, men's, closed toes, pair. S, M, L, XL. (add \$5.00)
- 8758 Convalescent boots, Hi-split style, ladies, open toes, pair. S M L
- 8759 Convalescent boots, Hi-split style, mens, open toes, pair. S M L XL (+\$5.00).
- 8761 Convalescent boots, Hi-split style, ladies', closed toes, pair, S, M, L,
- 8762 Convalescent boots, Hi-split style, men's, pair, S, M, L, XL. (add \$5.00)
- 8767 **Coccyx pads**, 6"x 6", pair.
- 8782 **Palm pads**, protect the palms from encroaching arthritic fingers, pair.
- 8784 **Palm cushion**, a double layer of sheepskin to keep the encroaching arthritic fingers away from the palm. Helps with sweat issues.
- 7130 **Woolskin detergent**, a specialty Tea Tree oil based wash for sheepskin & wool, with its anti-bacterial & anti-dust mite qualities. It has an additional ingredient to restore the suppleness to the leather. 500ml btl. with a built in measuring cup. Imported from Australia.
- 7135 Woolskin detergent, 5L jug. Imported from Australia.

Hi-temp “Green” medical sheepskin,

This medical grade sheepskin has had additional processing (compared to the regular grade medical sheepskin) to enable it to be washed and dried at higher temperatures. This includes being **washed** in the tannery to relax & restore the natural crimp, and improve the density of the fibres. It can be washed up to 80° C. This temperature makes it viable for **institutional** settings to maintain this product for their clients. It can be heat sanitized @ 72° enabling the product to be re-used with another client/patient. A prescribed washing chart will be provided.

Our regular medical grade can be laundered up to 40°C.

The natural characteristics of wool remain the same. Wool has the natural ability to absorb moisture up to 33% of its own weight. It cushions and provides comfort. Dense wool fibres assist in the distribution of weight and alleviate pressure points. Enjoy its benefits through **all** seasons. Air can circulate through the fibres, wicking the moisture to the atmosphere. Wool reduces shear (friction) problems. NB: Your initial washings of the Hi-temp sheepskins at high temperatures will see some nominal dye loss. Rinse in cold water afterwards.

Manufactured/distributed by:

THE WOOL TABLE SHEEPSKINS

Workshop: 80 Barbados Blvd, Unit 9, Scarborough, ON M1J 1K9

Mailing: P.O. Box 88036, Scarborough, ON M1M 3W1

416 755-9103 1 800 854-9774 Fax 647 351-5262

Please remember that The Wool Table Sheepskins is a specialty manufacturer of these premium quality products. This is all that we deal with. We can assist you with more technical details if needed. Our product lines are more extensive, as we try to attend to the variety of individual needs that this market demands. A major feature of our products is that each product has been designed to provide the most benefit for the user of it, e.g. we will not place a seam in a pressure sensitive area of a protective pad. Each one is thoroughly quality control checked. This minimizes the risk that one of your customers may bring back a defective product to you.

Graham McLeod, Owner.